



# Product Safety Recall



## Do you own this product?

### Scientrific — Vernier Go Direct Charging Station

The Vernier Go Direct Charging Station is designed to charge the Vernier Go Direct Wireless Educational Sensors.

Dates sold: 1 November 2017 – 15 December 2023



**Why the product is recalled:** A short circuit can occur in the charger when multiple wand-style GDX wireless educational sensors with damaged electrical contacts are plugged into the charging station. This can cause the charger to overheat and catch fire. Incidents have occurred.

**Hazard:** There is a risk of a serious burn injury and/or property damage if the charging station overheats and catches fire.



Consumers should stop using the product immediately and return them to Scientrific for repair, at no cost.



Consumers can contact Scientrific on

(02) 6645 8111 from 9:00 am to 5:00 pm, Monday to Friday,

or by email at: [info@scientrific.com.au](mailto:info@scientrific.com.au)



Purchased from: Scientrific Pty Ltd



Set yourself a reminder to check at home

- Take a photo of this notice
- Set a reminder on your phone

PRA 2024/20224

See [productsafety.gov.au](https://www.productsafety.gov.au) for  
Australian product recall information